

Advantages of GC-QQQ for the Analysis of Tobacco Alkaloids in Cigarette Filler

Joe Lisko

Research Chemist

jlisko@cdc.gov

ASMS Annual Conference


June 18, 2014

Tobacco Analysis Lab



Tobacco

- HPHCs (Harmful and Potentially Harmful Constituents)
– compounds that are toxic and/or carcinogenic
- Nicotine, TSNA's, PAH's, Metals, Flavors, Volatiles, Ammonia, TNCO



Smoke

- All Types of Products: Cigarettes, Cigars, Cigarillos, Roll-Your-Own, Smokeless Product, e-Cigarettes, etc.



Exposure

- Topography:
How people use the products

Tobacco and the Tobacco Matrix


- Tobacco use is the leading cause of preventable death worldwide – 5 million deaths annually



- The tobacco matrix is chemically complex and diverse.
- Approximately 4200 components have been found in tobacco – not including additives like flavors.
- Many factors affect the tobacco matrix – curing method, growing conditions, fertilizers, pesticides.
- Many classes of important compounds exist in the tobacco matrix: Metals, hydrocarbons, esters, lactones.....etc.

What are the minor tobacco alkaloids?

- Naturally occurring compounds in tobacco
- Nicotine-like in structure
- Make up 5% of total alkaloids in *N. tabacum*
 - nicotine accounts for 95%
- Biologically active in animals as well as humans
- Play a role in addiction and in carcinogenic TSNA formation




- Literature reported methods included GC-NPD, GC-FID, GC-MS and HPLC-MS

Method 1: GC-NPD

- Non-specific detection
- A matrix effect was observed despite using ITSD.
- Solid-Phase Microextraction (SPME) used for “pre-concentration” and injection
 - Partition coefficients changed with organic content
 - The SPME fiber gave variable results as it aged
 - Precision (RSD%) was 5-8%
 - Recovery was >94% for nicotine (other analytes not tested)

Alkaloids	Alkaloids/Internal Standard (peak ratio)	
	standards	spiked tobacco
Nornicotine	0.047	0.041
Myosmine	0.142	0.137
Anabasine	0.084	0.089


	Alkaloids/Internal Standard (peak ratio)			
	new fiber		aged fiber	
(Concentration, $\mu\text{g mL}^{-1}$)	10	20	10	20
Nornicotine	0.099	0.194	0.083	0.184
Myosmine	0.286	0.535	0.286	0.531
Anabasine	0.168	0.325	0.147	0.298



Method 2: GC-MS

- SPME sampling used (headspace)
- Used retention times and confirmation/quantitation ion ratios for more selective and accurate identification of analytes.
- Single Ion Monitoring (SIM) mode was used for quantitation.
- Initially attempted to use a denicotinized tobacco matrix for isotope dilution calibration.
 - Different relative response factors were observed between denicotinized tobacco and sample tobacco blends.
- Standard addition was used for calibration curves to minimize matrix effects between tobacco types
- Relatively high RSD and LOD's for the alkaloids.
- The method was applied to domestic and foreign cigarette filler.
 - Wide concentration ranges for all alkaloids were observed when comparing different international brands.
 - Most likely due to tobacco types used and/or quality of the tobacco.

Analyte	RSD (%)	Recovery (%)	LOD ($\mu\text{g/g}$)
nicotine	6.7	108	1.8
nor nicotine	9.9	100	0.134
anabasine	17.6	103	0.475
anatabine	20.6	108	0.209




GC-MS/MS Method

- Method required for rapid processing of large sample loads without loss of data integrity.
 - Through a collaboration with the FDA-CTP, CDC analyzes numerous tobacco products.
 - QC validation and protocol is always required.

Sample Preparation Improvements

- Tomato-leaf matrix used as a blank to reduce background.
 - Originally looked at Quest 3 denicotinized tobacco.
- Liquid-Liquid extraction
 - SPME inherently has more sample variability – due to more factors to overcome i.e. solubility and partition coefficients



GC-MS/MS Method - Continued

- **MS/MS Method Optimization – Agilent 7890 GC/7000 QQQ**
- **Isotopically labeled internal standards used for quantification (d_3 -nicotine and d_4 -nornicotine).**
- **Multiple Reaction Mode (MRM) used to quantify analytes**
 - Virtually eliminates background interferences
- **Validation (Precision/Accuracy) was done by spiking five samples (n=5) at three concentration levels.**
- **Improved precision, accuracy and lower LOD's than prior methods.**

Analyte	Precision, RSD (%)	Accuracy, Recovery (%)	LOD ($\mu\text{g/g}$)
nornicotine	0.9 – 3.3	102 – 112	0.08
myosmine	0.4 – 2.0	100 – 105	0.04
anabasine	0.8 – 1.7	98 – 107	0.12
anatabine	1.1 – 1.8	99 – 108	0.12
isonicotine	0.7 – 1.6	97 - 109	0.03

GC-MS/MS Method Application

- Analysis of cigarette filler from 50 popular US brands (n=7).

Analyte	Range ($\mu\text{g/g}$)	Mean ($\mu\text{g/g}$)	Median ($\mu\text{g/g}$)
Nornicotine	659 – 986	763	746
Myosmine	8.64 – 17.3	14.0	13.8
Anabasine	127 – 185	147	146
Anatabine	927 – 1390	1100	1090
Isonicotine	23.4 – 45.5	34.1	33.7

- Separated into mentholated and non-mentholated brands

Analyte	Non-Menthol (38 brands)			Menthol (12 brands)		
	Range ($\mu\text{g/g}$)	Mean ($\mu\text{g/g}$)	Median ($\mu\text{g/g}$)	Range ($\mu\text{g/g}$)	Mean ($\mu\text{g/g}$)	Median ($\mu\text{g/g}$)
Nornicotine	659 – 986	758	754	709 – 922	779	738
Myosmine	8.64 – 17.3	13.9	13.8	12.8 – 17.1	14.4	13.9
Anabasine	127 – 185	147	147	131 – 161	144	145
Anatabine	933 – 1390	1100	1090	927 – 1200	1100	1100
Isonicotine	23.4 – 45.5	33.8	33.7	30.0 – 45.5	35.0	33.4

Summary

Method	Number of Analytes	Sample Prep Method	QC (Y/N)	Validated (Y/N)	Internal Std (Y/N)	Precision (%)	Accuracy (%)	LOD (µg/g)
GC-NPD	4	SPME (in solution)	N	N	Y	5.0 – 8.0	>94*	ND
GC-MS	4	SPME (headspace)	Y	Y	Y#	6.7 – 20.6	100 - 108	0.134 – 1.8* (0.475)
GC-QQQ	5	Liquid-Liquid Extraction	Y	Y	Y#	0.4 – 3.3	97 - 112	0.03 – 0.12

ND = Not determined

* Nicotine

Isotopically labelled internal standards

Conclusions

- **GC-MS/MS approach provides a more efficient, more sensitive, and more robust method of alkaloids analysis**
 - LOD's were improved
 - Specificity was improved
 - Precision was improved
 - Sample throughput was increased
 - Tobacco matrix interferences were minimized





Acknowledgements

Stephen Stanfill – Team Lead, Smokeless Tobacco

Bryce Duncan

Cliff Watson – Lab Director, Tobacco Products

Ben Blount – Branch Chief, Tobacco Volatiles

FDA-Center for Tobacco Products - Funding

Thank You!

For more information please contact Centers for Disease Control and Prevention

1600 Clifton Road NE, Atlanta, GA 30333

Telephone: 1-800-CDC-INFO (232-4636)/TTY: 1-888-232-6348

Visit: www.cdc.gov | Contact CDC at: 1-800-CDC-INFO or www.cdc.gov/info

The findings and conclusions in this report are those of the authors and do not necessarily represent the official position of the Centers for Disease Control and Prevention.

National Center for Environmental Health
Division of Laboratory Science

